

Sentec Digital Monitor (SDM) Connectivity Overview

For a detailed overview of the SDM connectivity options including compatibility, technical information and instructions, refer to <https://www.sentec.com/transcutaneous-connectivity-overview/>

Patient Monitoring Systems (PMS)		
Ascom	Digistat® Smart Central	
Dräger	Infinity® Delta (XL), Infinity® Kappa, Infinity® Omega Serie	
GE Healthcare	Carescape B450/B650/B850, Solar 8000 M/i, Dash 5000	
Mindray	BeneView T5/T6/T8/T9, BeneVision N12, N15, N17, N19, N22	
Philips	IntelliVue MP40/50/60/70/80/90, MX400/450/500/550/600/700/750/800/850	
Spacelabs Healthcare	Flexport System Ultra View SL (90442A-03)	
Patient Data Management Systems (PDMS)		
Capsule	Medical Device Information Platform (MDIP)	
Dräger	Integrated Care Manager (ICM)	
iMDsoft	MetaVision	
Stowood	Visi-Download	
TZM	TZM (via) UMG to GE Healthcare's patient monitors	
Ventilators		
Air Liquide	Vendom 30/40	
Breas	Vivo 45 (LS/L), Nippy 4+	
EOVE	EO-150	
Philips	Trilogy Evo (O2), Trilogy EV300, A30 EFL, A40 Pro/EFL	
Thora Tech	Monsoon 4	
Digital Connectivity to Polygraphic and Polysomnographic Systems (PG/PSG)		
Cidelec	CID102L, CID102LB, CID-LXe, CID-LXe-206d, CID-LX, CID-LX-206d	
Nox Medical	NOX A1, NOX A1s (via C1), NOX C1, NOX T3	
Analog Connectivity to Polygraphic and Polysomnographic Systems (PG/PSG)		
Manufacturer/Brand	Device (**)	PSG Adapter Cable (Sentec Product Code)
Carefusion	SomnoStar Pro	PSG Cable E
Cidelec	CID102L8-D	PSG Cable P
Compumedics	E-Series, Grael (4k), Somte	PSG Cable E
	Siesta	PSG Cable G
Dr. Fenyves & Gut /Löwenstein Medical	SleepDoc Porti 7/8, MiniScreen PRO/plus	PSG Cable D
	Sonata	PSG Cable B
Natus	Comet / Aurora, Embla SDx, Embla Titanium (Isolation Box)	PSG Cable E
	Comlab, Deltamed Duo 44	PSG Cable H
	Deltamed Brain Box	PSG Cable J
	Embletta Gold, Embletta MPR PG	PSG Cable L
	Embla NDx, PAIB2 interface box	PSG Cable F
	Embla S7000 / N7000	PSG Cable A
	Embla Titanium	PSG Cable K
Mdrive, TX Proxy	PSG Cable N	
Neurovirtual	Sleepvirtual, BWII (SI1089)**, BWIII	PSG Cable E
SOMNOmedics	SOMNOscreen (plus) (SSC082)**	PSG Cable Q
Stowood Scientific Instruments	Visi-3, Visi-System	PSG Cable E
Philips Respironics	Alice 3 / 4 / 5 / 6 (square connector)	PSG Cable B
	Alice 3 / 4 / 5 / 6 (round connector)	PSG Cable F
Weinmann_MCC	SOMNOlab, SOMNOlab2	PSG Cable C

** Product number of accessory cables or optocoupler required for installation in brackets behind the product name.

Note: This is a global overview and may contain information on products not released in all markets/countries. Information is subject to change without notice.